

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Jennifer Shirley - AZDPS AGENCY LHC PD

DATE 10/12/11 TIME 1012

INTOXILYZER SERIAL # 80-001029 LOCATION LHC PD Jail

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Jennifer Shirley

INTOXILYZER 8000

Location: LHCPD MAIN JAIL

Serial Number: 80-001029

Core Version: 8105.48

10/12/2011

10:16:46

Standard Lot#: 35710100A1

Last Changed By: K. DRELLER #121

QAS: J. SHIRLEY #7182

AZ DPS

Operator: J. SHIRLEY #7182

AZ DPS

Test	g/210L	Time
Air Blank	0.000	10:17:32
0.100 Cal Check	0.099	10:17:53
Air Blank	0.000	10:18:22

Preventative Maintenance
Performed -- Timer Reset